

# Sample Environment for Soft Matter at ILL SANS Instruments

Lionel PORCAR, Ralf SCHWEINS, Sylvain PREVOST, Olga MATSARSKAIA, Anne MARTEL

Institut Laue Langevin, Grenoble, France

Various techniques can be combined *in situ* with SANS either to condition the sample or to obtain complementary information. This talk will give an overview of the ones which are available at ILL SANS instruments, as well as some set up developed by our users.

Among sample conditioners, dialysis cell, rheometer, high-pressure cell, chromatography and last developments will be presented.

As complementary probes, SAXS, DLS and fluorimetry will be highlighted.